

H3N2 HA (His Tag) recombinant protein



Catalog Number: 502542

General Information

Gene Name Synonym

Hemagglutinin HA1 chain; Hemagglutinin HA2 chain

Protein Construction

A DNA sequence encoding the Influenza A virus (A/Aichi/2/1968 (H3N2)) hemagglutinin (P03437.1) (Met1-Trp530), termed as HA, was expressed with a polyhistidine tag at the C-terminus.

Organism

H3N2

Expression Host

Baculovirus-Insect Cells

QC Testing

Activity

Measured by its ability to agglutinate guinea pig red blood cells.

HA titer is 0.1-0.5µg/mL for 1%GRBC

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 17

Molecular Mass

The recombinant hemagglutinin of Influenza A virus (A/Aichi/2/1968 (H3N2)) consists 525 amino acids and predicts a molecular mass of 59 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10 % glycerol, pH 8.0.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Influenza A H3N2 (A/Aichi/2/1968) Hemagglutinin / HA Protein (His Tag)