

# Vitamin D binding protein antibody



Catalog Number: 107542

---

## Product name

Vitamin D binding protein antibody

## Specificity

Human; other species not tested.

## Antibody description

Vitamin D binding protein Mouse Monoclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 52-58 kDa

## Preparation

This antibody was obtained by immunization of Vitamin D binding protein recombinant protein (Accession Number: NM\_000583). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107542(VDBP,GC Antibody) at dilution of 1:4000