

AKR1C4 antibody



Catalog Number: 107942

Product name

AKR1C4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

AKR1C4 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, HepG2 cells, mouse liver tissue. Positive IF detected in HepG2. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of AKR1C4 recombinant protein (Accession Number: NM_001818). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107942(AKR1C4 antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using AKR1C4 antibody Catalog No:107942 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).