

# Anti-FOXA1 Rabbit antibody



Catalog Number: 168942

## Product name

Anti-FOXA1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FOXA1 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human FOXA1

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions


WB: 1/1000-1/5000

IHC: 1/20

ICC/IF: 1/20

FC: 1/20

## Validations



Western blot detection of FOXA1 in C6 cell lysates using FOXA1 Rabbit mAb(1:1000 diluted).Predicted band size:49kDa.Observed band size:49kDa.