

# DIMT1L antibody



Catalog Number: 109942

## Product name

DIMT1L antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DIMT1L Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, Jurkat cells, Raji cells. Observed molecular weight by Western-blot: 35kd

## Preparation

This antibody was obtained by immunization of DIMT1L recombinant protein (Accession Number: NM\_001348076). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:109942(DIMT1L antibody) at dilution of 1:500