

Anti-CD58 antibody



Catalog Number: 141242

Product name

Anti-CD58 antibody

Immunogen

[Human LFA-3 / CD58 \(His Tag\) recombinant protein](#)

Specificity

Human LFA-3 / CD58

Antibody description

Mouse Monoclonal to Human CD58

Preparation

This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human LFA-3 / CD58 and conjugated with PerCP under optimum conditions, the unreacted PerCP was removed.

Formulation

Aqueous solution containing 0.5% BSA and 0.09% sodium azide

Storage

This antibody is stable for 12 months from date of receipt when stored at 2°C-8°C. Protected from

prolonged exposure to light. Do not freeze ! Sodium azide is toxic to cells and should be disposed of properly. Flush with large volumes of water during disposal.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

FCM

Dilutions

FCM

Validations



Human LFA-3 / CD58 Flow Cytometry (FC) 15900