

GTPBP2 antibody

Catalog Number: 111242

Product name

GTPBP2 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human,Mouse,Rat; other species not tested.

Antibody description

GTPBP2 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human spleen tissue. Positive IP detected in mouse spleen tissue. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of GTPBP2 recombinant protein (Accession Number: BC028347). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:111242(GTPBP2 antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111242(GTPBP2 antibody) at dilution of 1:50



(under 10x lens)

IP Result of anti-GTPBP2 (IP:Catalog No:111242, 4ug; Detection:Catalog No:111242 1:300) with

GTPBP2 antibody



Catalog Number: 111242

mouse spleen tissue lysate 4000ug.