

# PRR14 antibody



Catalog Number: 114242

## Product name

PRR14 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PRR14 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HL-60 cells, mouse brain tissue. Observed molecular weight by Western-blot: 64 kDa

## Preparation

This antibody was obtained by immunization of PRR14 recombinant protein (Accession Number: NM\_024031). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114242(PRR14 Antibody) at dilution of 1:1000