

# Anti-XRCC1 (3F6) Mouse antibody



Catalog Number: 168242

## Product name

Anti-XRCC1 (3F6) Mouse antibody

## Ig Type

Mouse IgG1

## Specificity

Human

## Applications

WB

## Antibody description

XRCC1 (3F6) Mouse monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: Purified recombinant human XRCC1 protein fragments expressed in E.coli.

## Validations

## Formulation


Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal



Western blot detection of XRCC1 in Raji, K562, HeLa, Jurkat and THP-1 cell lysates using XRCC1 mouse mAb (1:1000 diluted). Predicted band size: 82KDa. Observed band size: 82KDa.