

CEA antibody

Catalog Number: 109242

Product name

CEA antibody

Specificity

Human; other species not tested.

Antibody description

CEA Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human colon tissue, human lung cancer tissue, human pancreas cancer tissue, human stomach cancer tissue.

Preparation

This antibody was obtained by immunization of CEA recombinant protein (Accession Number: NM_004363). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109242(CEA antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109242(CEA antibody) at dilution of 1:50 (under 40x lens)