

Human UBE2T/HSPC150/PIG50 (His tag) recombinant protein



Catalog Number: 519242

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 T, Cell Proliferation-Inducing Gene 50 Protein, Ubiquitin Carrier Protein T, Ubiquitin-Protein Ligase T, UBE2T

Protein Construction

A DNA sequence encoding the Human UBE2T/HSPC150/PIG50 Met1-Val197 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2T/HSPC150/PIG50 (His tag)
recombinant protein