

IκB beta antibody



Catalog Number: 111642

Product name

IκB beta antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IκB beta Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue, human endometrial cancer tissue. Positive WB detected in HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of IκB beta recombinant protein (Accession Number: NM_002503). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111642(NFKBIB antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111642(NFKBIB antibody) at dilution of 1:100 (under 10x lens)