

MEK1-Specific antibody



Catalog Number: 112642

Product name

MEK1-Specific antibody

Specificity

Human; other species not tested.

Antibody description

MEK1-Specific Rabbit Polyclonal antibody.
Positive IHC detected in human endometrial cancer tissue, human colon tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002755).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using Catalog No:112642(MEK1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue slide using Catalog No:112642(MEK1 Antibody) at dilution of 1:50 (under 40x lens)