

# PDXK antibody

Catalog Number: 113642

## Product name

PDXK antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PDXK Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in human liver tissue, HeLa cells. Observed molecular weight by Western-blot: 40kd, 35kd

## Preparation

This antibody was obtained by immunization of PDXK recombinant protein (Accession Number: NM\_003681). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113642(PDXK antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113642(PDXK antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog

# PDXK antibody



Catalog Number: 113642

---

No:113642(PDXK antibody) at dilution of 1:50 (under 40x lens)