

# DIO1 Polyclonal Antibody



Catalog Number: 166642

## Product name

Rabbit IgG

DIO1 Polyclonal Antibody

## Applications

## Specificity

WB

Human

## Dilutions

## Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to DIO1

## Validations

## Preparation

Antigen: A synthetic peptide of human DIO1

## Formulation

PBS with 0.02% sodium azide, pH7.3.


## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type



Western blot - DIO1 Polyclonal Antibody

Western blot analysis of extracts of HepG2 cells, using DIO1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.