

# Anti-ALB antibody



Catalog Number: 176642

## Product name

Anti-ALB antibody

## Specificity

Bovine

## Antibody description

Mouse monoclonal antibody to ALB

## Preparation

This antigen of this antibody was protein

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB

## Dilutions

WB: 1:500-1:1,000

## Validations




Fig1: Western blot analysis on BSA using anti-BSA Mouse mAb (Cat. # 176642#).