

Anti-ALB antibody

Catalog Number: 176642

Product name

Anti-ALB antibody

Specificity

Bovine

Antibody description

Mouse monoclonal antibody to ALB

Preparation

This antigen of this antibody was protein

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB

Dilutions

WB: 1:500-1:1,000


Validations

Fig1: Western blot analysis on BSA using anti-BSA Mouse mAb (Cat. # 176642#).