

Anti-Rb Rabbit antibody



Catalog Number: 168642

Product name

Anti-Rb Rabbit antibody

Specificity

Human, Mouse

Antibody description

Rb Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of mouse Rb

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of Rb in C6, Jurkat, HeLa cell lysates using Rb Rabbit mAb(1:1000 diluted). Predicted band size: 106kDa. Observed band size: 106kDa.