

CXCR1 antibody

Catalog Number: 109642

Product name

CXCR1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CXCR1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, mouse placenta tissue. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000634). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109642(CXCR1 antibody) at dilution of 1:300



IP Result of anti-CXCR1 (IP:Catalog No:109642, 4ug; Detection:Catalog No:109642 1:500) with Jurkat cells lysate 4000ug.