

# Cynomolgus IL-21R (His Tag) recombinant protein



Catalog Number: 500742

## General Information

### Protein Construction

A DNA sequence encoding the cynomolgus IL21R (XP\_005591578.1) (Met1-Pro236) was expressed with a polyhistidine tag at the C-terminus.

### Organism

Cynomolgus

### Expression Host

Human Cells

## QC Testing

### Purity

> 90 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Cys 20

### Molecular Mass

The recombinant cynomolgus IL21R comprises 228 amino acids and has a calculated molecular mass of 26.4 KDa. The apparent molecular mass of

it is approximately 44.5 KDa in SDS-PAGE.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide

### Storage


Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ .

It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Cynomolgus IL-21R / Interleukin-21 Receptor Protein (His Tag) SDS-PAGE