

# Human GAS7/KIAA0394 (His tag) recombinant protein



Catalog Number: 518742

## General Information

### Gene Name Synonym

Growth Arrest-Specific Protein 7, GAS-7, GAS7, KIAA0394

### Protein Construction

A DNA sequence encoding the Human GAS7/KIAA0394 Met1-Ile412 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 10% Glycerol, pH 8.8.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

## SDS-PAGE



Human GAS7/KIAA0394 (His tag) recombinant protein