

NTF3 Polyclonal Antibody



Catalog Number: 164052

Product name

NTF3 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to NTF3

Preparation

Antigen: Recombinant protein of human NTF3

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000


IHC 1:100 - 1:200

Validations




Western blot - NTF3 Polyclonal Antibody

Western blot analysis of extracts of mouse heart, using NTF3 antibody at 1:250 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - NTF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using NTF3 antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - NTF3 Polyclonal Antibody

NTF3 Polyclonal Antibody



Catalog Number: 164052

Immunohistochemistry of paraffin-embedded

mouse kidney using NTF3 antibody at dilution of 1:100 (40x lens).