

# POLR3G antibody

Catalog Number: 114052

## Product name

POLR3G antibody

## Specificity

Human; other species not tested.

## Antibody description

POLR3G Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells, HEK-293 cells, K-562 cells. Observed molecular weight by Western-blot: 26-30 kDa

## Preparation

This antibody was obtained by immunization of POLR3G recombinant protein (Accession Number: XM\_011543100). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114052(POLR3G Antibody) at dilution of 1:1000



IP Result of anti-POLR3G (IP:Catalog No:114052, 3ug; Detection:Catalog No:114052 1:500) with HeLa cells lysate 2800ug.