

Anti-LBP antibody



Catalog Number: 176052

Product name

Anti-LBP antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to LBP

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human lbp. aa range:221-270

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using Antibody diluted at 1:1000. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).