

CD206 antibody



Catalog Number: 107052

Product name

CD206 antibody

Specificity

Human; other species not tested.

Antibody description

CD206 Mouse Monoclonal antibody. Positive WB detected in human liver tissue, human placenta tissue. Observed molecular weight by Western-blot: 170 kDa

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_002438). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2a

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:107052(MRC1 Antibody) at dilution of 1:1000