

Anti-CaM Kinase 4/CaMKIV antibody



Catalog Number: 100452

Product name

Anti-CaM Kinase 4/CaMKIV antibody

Immunogen

[Mouse CaM Kinase 4/CaMKIV \(His & GST Tag\) recombinant protein](#)

Specificity

Mouse CaM Kinase 4 / CaMKIV

Antibody description

Rabbit monoclonal to CaM Kinase 4/CaMKIV

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Mouse CaM Kinase 4 / CaMKIV (rM CaM Kinase 4 / CaMKIV; P08414; Met1-Tyr469).

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

IHC-P

Dilutions


IHC-P: 1-10µg/mL

Validations



CaM Kinase 4 / CaMKIV Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse CAMK4 in mouse prostate with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded



sections).

CaM Kinase 4 / CaMKIV Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of mouse CAMK4 in mouse brain with rabbit monoclonal antibody (2 µg/mL, formalin-fixed paraffin embedded sections).