

FAF2 antibody

Catalog Number: 110452

Product name

FAF2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FAF2 Rabbit Polyclonal antibody. Positive IF detected in PC-3 cells. Positive IHC detected in human testis tissue. Positive WB detected in mouse brain tissue, HeLa cells, HepG2 cells, human brain tissue, human heart tissue, human pancreas tissue, mouse heart tissue, mouse lung tissue, mouse pancreas tissue, mouse testis tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 53 kDa

Preparation

This antibody was obtained by immunization of FAF2 recombinant protein (Accession Number: NM_014613). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations





mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110452(FAF2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human testis using Catalog No:110452(FAF2 antibody) at dilution of 1:50 (under 10x lens)

FAF2 antibody

Catalog Number: 110452



IP Result of anti-FAF2 (IP:Catalog No:110452, 3ug; Detection:Catalog No:110452 1:1000) with mouse brain tissue lysate 5000ug.

Immunofluorescent analysis of PC-3 cells using Catalog No:110452(FAF2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)