

HPS3 antibody

Catalog Number: 111452

Product name

HPS3 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

HPS3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HEK-293T cells. Positive IHC detected in human colon cancer tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 113 kDa



Preparation

This antibody was obtained by immunization of HPS3 recombinant protein (Accession Number: BC022062). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111452(HPS3 antibody) at dilution of 1:300

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111452(HPS3 antibody) at dilution of 1:100 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:1000



Immunofluorescent analysis of HeLa cells using

HPS3 antibody



Catalog Number: 111452

Catalog No:111452(HPS3 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-rabbit IgG(H+L)