

MAGEF1 antibody

Catalog Number: 112452

Product name

MAGEF1 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEF1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of MAGEF1 recombinant protein (Accession Number: NM_022149). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112452 (MAGEF1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112452(MAGEF1 Antibody) at dilution of 1:162 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:112452(MAGEF1 Antibody) at dilution of 1:162 (under 40x lens)