

# PAX1-Specific antibody



Catalog Number: 107452

---

## Product name

PAX1-Specific antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PAX1-Specific Mouse Monoclonal antibody.  
Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006192).  
Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:107452(PAX1-Specific Antibody) at dilution of 1:1000