

BLNK antibody

Catalog Number: 108452

Product name

BLNK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

BLNK Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 68 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of BLNK recombinant protein (Accession Number: NM_013314). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:117187(BLNK antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:117187(BLNK antibody) at dilution of 1:100 (under 10x lens)