

Anti-ANKRD30A antibody



Catalog Number: 175852

Product name

Anti-ANKRD30A antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ANKRD30A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human an30a. aa range:1041-1090

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using ANKRD30A Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).