

PPP1R14D antibody



Catalog Number: 114152

Product name

PPP1R14D antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PPP1R14D Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of PPP1R14D recombinant protein (Accession Number: NM_017726). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:114152(PPP1R14D Antibody) at dilution of 1:1000