

Adecatumumab reference antibody



Catalog Number: 177152

Product name

Monoclonal

Adecatumumab reference antibody

Ig Type

Specificity

IgG1

Human EpCAM / TROP1 / CD326

Applications

Antibody description

ELISA

Anti-EpCAM / TROP1 / CD326 Antibody
(Adecatumumab)

Validations

Preparation

Recombinant expression and purified from CHO cells.


Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality



SDS-PAGE

Anti-EpCAM / TROP1 / CD326 Antibody (adecatumumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.