

UFSP2 antibody

Catalog Number: 116552

Product name

UFSP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UFSP2 Rabbit Polyclonal antibody. Positive IF detected in HeLa cells. Positive IHC detected in human liver tissue, human brain tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human spleen tissue, human testis tissue. Positive IP detected in HepG2 cells. Positive WB detected in mouse pancreas tissue, HeLa cells, HepG2 cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 53kd

Preparation

This antibody was obtained by immunization of UFSP2 recombinant protein (Accession Number: NM_018359). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP, IF

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:10-1:100

Validations




mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:116552(UFSP2 antibody) at dilution of 1:400




Immunohistochemical of paraffin-embedded human liver using Catalog No:116552(UFSP2 antibody) at dilution of 1:100 (under 10x lens)

UFSP2 antibody


Catalog Number: 116552



Immunohistochemical of paraffin-embedded human liver using Catalog No:116552(UFSP2 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-UFSP2 (IP:Catalog No:116552, 3ug; Detection:Catalog No:116552 1:1000) with HepG2 cells lysate 1700ug.



Immunofluorescent analysis of HeLa cells using Catalog No:116552(UFSP2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)