

Anti-FOXD3 (3D11) Mouse antibody



Catalog Number: 168552

Product name

Anti-FOXD3 (3D11) Mouse antibody

Specificity

Human, Mouse, Monkey

Antibody description

FOXD3 (3D11) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human FOXD3 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications

WB, IHC-P


Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ELISA: 1/10000

Validations



Immunohistochemical analysis of paraffin-embedded lung cancer tissues (left) and ovarian cancer tissues (right) using FOXD3 mouse mAb with DAB staining.



with DAB staining.

Western blot detection of FoxD3 in MCF7 and COS7 cell lysates using FoxD3 mouse mAb (1:1000 diluted). Predicted band size: 48KDa. Observed band size: 48KDa.