

# Anti-FOXD3 (3D11) Mouse antibody



Catalog Number: 168552

## Product name

Anti-FOXD3 (3D11) Mouse antibody

## Specificity

Human, Mouse, Monkey

## Antibody description

FOXD3 (3D11) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human FOXD3 expressed in E. Coli.

## Formulation

Ascitic fluid containing 0.03% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications

WB, IHC-P


## Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000

ELISA: 1/10000

## Validations



Immunohistochemical analysis of paraffin-embedded lung cancer tissues (left) and ovarian cancer tissues (right) using FOXD3 mouse mAb



with DAB staining.

Western blot detection of FoxD3 in MCF7 and COS7 cell lysates using FoxD3 mouse mAb (1:1000 diluted). Predicted band size: 48 kDa. Observed band size: 48 kDa.