

# NARFL-Specific antibody

Catalog Number: 112952

## Product name

NARFL-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NARFL-Specific Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IHC detected in human kidney tissue, human heart tissue. Observed molecular weight by Western-blot: 26-28 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_022493). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations




HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112952(NARFL-Specific antibody) at dilution



of 1:800

Immunohistochemical of paraffin-embedded human kidney using Catalog No:112952(NARFL-Specific antibody) at dilution of 1:100 (under 40x



lens)

Immunohistochemical of paraffin-embedded human kidney using Catalog No:112952(NARFL-Specific antibody) at dilution of 1:100 (under 10x lens)