

Anti-GATA1 Rabbit antibody



Catalog Number: 168952

Product name

Anti-GATA1 Rabbit antibody

Specificity

Human

Antibody description

GATA1 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human GATA1

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, ChIP

Dilutions

WB: 1/2000


IHC: 1/20

ICC/IF: 1/50

IP: 1/20

ChIP: 1/20

Validations



Western blot detection of GATA1 in K562, HeLa cell lysates using GATA1 Rabbit mAb(1:1000 diluted). Predicted band size: 43kDa. Observed band size: 43kDa.