

# DNAJC7 antibody

Catalog Number: 109952

## Product name

DNAJC7 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

## Specificity


Human, Mouse; other species not tested.

IF: 1:10-1:100

## Antibody description

DNAJC7 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, mouse testis tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human pancreas cancer tissue. Positive FC detected in HepG2 cells. Observed molecular weight by Western-blot: 56 kDa

## Validations




## Preparation

This antibody was obtained by immunization of DNAJC7 recombinant protein (Accession Number: NM\_003315). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109952(DNAJC7 antibody) at dilution of 1:500

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109952(DNAJC7 antibody) at dilution of 1:50 (under 10x lens)

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IHC, IF, FC

## Dilutions

Recommended Dilution:




# DNAJC7 antibody



Catalog Number: 109952

Immunofluorescent analysis of HepG2 cells, using DNAJC7 antibody Catalog No:109952 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



1X10<sup>6</sup> HepG2 cells were stained with 0.2ug DNAJC7 antibody (Catalog No:109952, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.