

NMNAT1 antibody



Catalog Number: 113252

Product name

NMNAT1 antibody

Specificity

human; other species not tested.

Antibody description

NMNAT1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in human testis tissue, human brain tissue, human heart tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of NMNAT1 recombinant protein (Accession Number: XM_047428076). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113252(NMNAT1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:113252(NMNAT1 antibody) at dilution of 1:50 (under 10x lens)