

Human AOC3/VAP1 (His tag) recombinant protein



Catalog Number: 518252

General Information

Gene Name Synonym

Membrane primary amine oxidase; Copper amine oxidase; HPAO; Semicarbazide-sensitive amine oxidase; SSAO; Vascular adhesion protein 1; VAP-1; AOC3; VAP1

Protein Construction

A DNA sequence encoding the Human AOC3/VAP1 Arg27-Asn763 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 500mM NaCl, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human AOC3/VAP1 (His tag) recombinant protein