

# IL 1 alpha antibody



Catalog Number: 111652

## Product name

IL 1 alpha antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL 1 alpha Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human placenta tissue, Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 31 kDa, 17 kDa

## Preparation

This antibody was obtained by immunization of IL 1 alpha recombinant protein (Accession Number: NM\_000575). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111652(IL1A Antibody) at dilution of 1:600