

# ANG Polyclonal Antibody



Catalog Number: 162652

---

## Product name

ANG Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ANG

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-147 of human ANG (NP\_001091046.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



Immunofluorescence - ANG Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ANG antibody . Blue: DAPI for nuclear staining.