

ANG Polyclonal Antibody



Catalog Number: 162652

Product name

ANG Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ANG

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 25-147 of human ANG (NP_001091046.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Immunofluorescence - ANG Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ANG antibody . Blue: DAPI for nuclear staining.