

# PDCL antibody



Catalog Number: 113652

## Product name

PDCL antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

## Immunogen

[Human PDCL Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PDCL Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 40-42 kDa

## Preparation

This antibody was obtained by immunization of PDCL recombinant protein (Accession Number: NM\_005388). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113652(PDCL antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human placenta using Catalog No:113652(PDCL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:113652(PDCL antibody) at dilution of 1:50 (under 10x lens)

# PDCL antibody



Catalog Number: 113652

---

antibody) at dilution of 1:50 (under 40x lens)