

# RHAMM,Lgals9 antibody



Catalog Number: 114652

## Product name

RHAMM,Lgals9 antibody

## Specificity

mouse, rat; other species not tested.

## Antibody description

RHAMM,Lgals9 Rabbit Polyclonal antibody.  
Positive WB detected in NIH/3T3 cells. Observed  
molecular weight by Western-blot: 32kd,70kd

## Preparation

This antibody was obtained by immunization of  
Peptide (Accession Number: NM\_012964).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



NIH/3T3 cells were subjected to SDS PAGE  
followed by western blot with Catalog  
No:114652(RHAMM,Lgals9 antibody) at dilution  
of 1:400