

# F5 Polyclonal Antibody



Catalog Number: 161352

## Product name

F5 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to F5

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 300-510 of human F5 (NP\_000121.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



Western blot - F5 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using F5 antibody .