

# SLCO4A1 antibody

Catalog Number: 115352

## Product name

SLCO4A1 antibody

## Specificity

Human; other species not tested.

## Antibody description

SLCO4A1 Rabbit Polyclonal antibody. Positive IP detected in L02 cells. Observed molecular weight by Western-blot: 77 kDa

## Preparation

This antibody was obtained by immunization of SLCO4A1 recombinant protein (Accession Number: NM\_016354). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP

## Dilutions

Recommended Dilution:

IP: 1:200-1:2000

## Validations



IP Result of anti-SLCO4A1 (IP:Catalog No:115352, 4ug; Detection:Catalog No:115352 1:500) with L02 cells lysate 1800ug.