

Anti-TMC8 antibody



Catalog Number: 176352

Product name

Anti-TMC8 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TMC8

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tmc8. aa range:601-650

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of HeLa cells using TMC8 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000